Nit-Occlud® PDA-R

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Instructions for Use

Nit-Occlud® PDA-R

Device for the closure of a patent ductus arteriosus

Product description

The device is built from one nitinol wire without soldering (Fig. 01 - 03). The device has a unique configuration to fix the implant to the aortic ampulla. The cylindrical structure expands against the wall of the duct without protrusion into the Aorta or the Pulmonary arteries.


Fig. 01: Distal disc view

Fig. 02: Lateral view


Fig. 03: Proximal disc view

The Nit-Occlud® PDA-R utilizes its elasticity to adapt to a wide range of different ductal morphologies. Several polyester membranes are fixed inside the implant to promote thrombosis thus accelerating closure. Also, the implant has two platinum markers for better radiological visibility.

The device is attached to a controlled release system (Fig. 4 and 5, Table 1). Nit-Occlud® Implantation Sheaths preferable with an internal diameter between 6 to 9 French are used for the deployment of the device.

The Nit-Occlud® PDA-R is a Class III cardiac medical device (according to Directive 93/42/ECC, Annex IX, Rule 8). It is implanted using a minimally invasive technique called cardiac catheterization.


Fig. 04: Lateral view


Fig. 05: Complete system

Letter	Description
D	Device
P	Pusher*
A	Adaptor
TS	Transportation Sheath
PC	Pusher Catheter
Y	"Y" connector
SS	Security seal
T	Torque device
W	Lock wire
d	Dilator

Tbl. 01

*The bigger versions of the Nit-Occlud® PDA-R (Ref. No.160106, 160107, 160108) are delivered with a pusher catheter, which improves stability. The pusher catheter is of light blue color.

General Information

The implantation of the devices with cardiac catheterization, must be exclusively carried out by a trained physician. It is recommended to read thoroughly these instructions and consider the detailed characteristics on the label.

Product Identification

Each product label has peel-off labels to allow the product to be identified precisely. Preserve the information of the product; remove the adhesive label and stick it to the medical records of each patient.

MR-Compatibility

MRI labeling based on the test results:

The Nit-Occlud® PDA-R was determined to be MR-conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-08. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania.

Non-clinical testing demonstrated that the Nit-Occlud® PDA-R Occluder is MR Conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3 Tesla or less
- Maximum spatial gradient magnetic field of 720 Gauss/cm or less

MRI-Related Heating

In non-clinical testing, the Nit-Occlud® PDA-R Occluder produced the following temperature rise during MRI performed for 15 min of scanning (i.e., per pulse sequence) in the 3 Tesla (3 Tesla/128 MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI)

MR system: Highest temperature change +1.7°C

Therefore, the MRI-related heating experiments for the Nit-Occlud® PDA-R Occluder at 3 Tesla using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 2.9 W/kg (i.e., associated with a calorimetry measured whole body averaged value of 2.7 W/kg)

indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.7°C.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Nit-Occlud® PDA-R. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

Pulse Sequence	T1 - SE	T1 - SE	GRE	GRE
Signal Void Size	145 mm ²	196 mm ²	534 mm ²	448 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Indication

The Nit-Occlud® PDA-R is an implant developed for the transcatheter closure of the Patent Ductus Arteriosus (PDA), in defects between 2 and 8 mm (pulmonary orifice).

Contraindications

Pathological or physical conditions precluding the implantation of a Nit-Occlud® PDA-R, such as:

- Pregnancy
- Patients with a weight lower than 10 kg
- Associated cardiomyopathy
- Immunologic or oncologic diseases
- Active infections
- Febrile syndrome with unknown origin
- Pulmonary hypertension higher than 2/3 of the systemic pressure
- Hematological diseases
- Allergy to the device's components (nickel)
- Ductus with atypical morphology, calcified or aneurysmal
- Endocarditis
- Patient with the Eisenmenger syndrome
- Allergy to the contrast used for the catheterization
- Existing thrombosis anywhere along the delivery catheter course
- Patients whose veins are too small or have anomalous connections not allowing safe or direct access to the right heart

Possible Complications

- Embolization of Nit-Occlud® PDA-R
- Obstruction of the aorta and/or pulmonary arteries
- Hemolysis
- Air embolism stroke
- Thrombosis stroke
- Weak or absent peripheral pulse
- Vascular perforation
- Hemorrhage
- Hematoma in the puncture place
- Arrhythmia
- Fever
- Endocarditis
- Death

Precautions

WARNINGS

-  **Perform an angiography to choose the proper device**
-  **The patient has to be heparinized in order to obtain an activated clotting time (ACT value) of at least 200 to 250 seconds at the time of the implantation.**
-  **Do not inject contrast media through Nit-Occlud® PDA-R System.**
-  **If the device doesn't configure correctly or protrudes either into the aorta or pulmonary artery, retrieve it and repeat the procedure with an adequate device.**

- Only use the Nit-Occlud® PDA-R for occlusion of a patent ductus. Do not use this device for other pathologies.
- Verify that the device corresponds to the angiographic measurements of the defect. Ensure that the Implantation Sheath is suitable for the chosen device. Follow these steps identically on the instrumentation table before the procedure.
- In case of requiring a very small device in a patient over 40 kg, it is recommended to use a stabilizing sheath, due to the contraction of the RV which makes the catheter of the implant unstable.
- Ensure that the Implantation Sheath is appropriate for the patient's weight to prevent damage to the vessels.
- To avoid the risk of air embolization, all catheters must be carefully flushed.
- During the procedure keep the system straight. Avoid to pull the lock wire accidentally: Do not manipulate the seal from the Y connector or the torque device until you are sure that it is the right moment to release the device.
- Do not manipulate the retaining wires which are fixed to the Y connector.
- Avoid exaggerated manipulation of the device, extremely long procedures or insufficient heparinization of the patient.
- Should the implant embolize, attempt interventional capturing. Inform yourself about the capturing techniques and procure the necessary equipment (e.g. snares, forceps, long sheath).
- Clinical evaluation indicates that in 0.09% of the patients could develop complications causing the need for surgery.
- The literature informs that the mortality is lower than one per one thousand cases with similar devices.
- Each device has been packaged under a controlled environment and sterilized with ethylene oxide. It is provided in a pyrogen-free condition. The device is for single use only. It is recommended that the product is stored in a dry place, protected against extreme humidity. The products should not be exposed to direct sun light.
- Do not use the product if the package is damaged, opened, wet and/or if the expiration date has passed or if some components of the system are damaged. In any of these cases please contact the manufacturer or distributor.
- Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization of single-use devices may result in degraded performance or a loss of functionality. Reuse of single-use devices may result in exposure to pathogens such as viruses, bacteria, fungi, or prions.
- The residuals of medical products and their components can represent a biological danger. Hence, the norms established in each region must be followed to dispose medical waste.

Patient's Medication

As antibacterial prophylaxis administer a wide spectrum antibiotic, one dose before and one to two doses after implantation is suggested or according to your institution.

It is essential to provide the patient with an efficient endocarditis prophylaxis until six months after implantation.

Ductus Measurement Parameters and Choosing the Device

Make an aortogram in lateral projection:

- A Measure the minimal diameter (md) of the ductus.
- B Measure the aortic ampulla diameter (am).
- C Measure the length (ld) between the pulmonary and aortic orifice. Fig. 6.


Fig. 06


Fig. 07


Fig. 08

Letter	Description
PA	Pulmonary artery
Ao	Aorta
md	Minimum diameter
am	Ampulla diameter
ld	Length of ductus

Tbl. 02

Choosing the device

$S = md \times (1,5 \text{ to } 2)$

The selection of the device is dictated by the minimum diameter of the defect:

The diameter (S) of the device should be 2 times bigger than the minimum diameter (md) of the defect. It must never and in no case be smaller than 1.5 times of the defect's diameter. The physician must decide the size of the device based on morphology of the defect and patient's weight according to the suggestions in Tab. 03.

The selected device should not cause obstruction in the Aorta or Pulmonary artery - with hemodynamic relevance - due to its length (L) or the diameter of the retention disk (D). Fig. 6 - 8.

Use the table to determine the recommended Implantation Sheath and the catalog number to order the product.

Reference number	Suggested prescription for PDA Diameter (mm)	Stent (mm) S	Disc (mm) D	Length (mm) L	Nit-Occlud® Implantation Sheath (F)
160102	2	4	8	6,5	6
160103	3	5,5	10	7	6
160104	4	7	12	8,5	6
160105	5	8,5	14	9,5	6
160106	6	10	16	11	8
160107	7	11,5	18	12	8
160108	8	13	20	13,5	9

Tbl. 03

The sizes/diameters above relate to the use of the Nit-Occlud® implantation sheath.

If using other sheaths, the size/diameter requirements may differ.

Description of the Implantation Technique

- The recommended Nit-Occlud® Implantation Sheath should be positioned in standard fashion across the ductus from the venous route.
- Verify visually the completeness of the device and check if all components are in place.
- Pull the pusher (P), until the device (D) is completely inside the transportation sheath (TS) (Fig. 09).


Fig. 09

- Flush the whole system with heparinised solution (11.U./1ml), according to the usual technique. It is recommended to aspirate the heparinized solution at least once to eliminate air bubbles.
- Screw the transportation sheath (TS) to the Implantation Sheath (IS) via the luer connector (Fig. 10).


Fig. 10

- In case you use an Implantation Sheath (IS) (max. Length = 90cm), with a hemostatic valve, introduce the adaptor (A) into the Implantation Sheath (IS), and remove the dilator (d) (Fig. 11-12).
- Connect the transportation sheath (TS) to the adaptor (A) (Fig. 13).


Fig. 11


Fig. 12


Fig. 13

- Advance the pusher (P) until the device reaches the end of the Implantation Sheath (IS). Fig. 14.


Fig. 14


Fig. 15

- Verify the position of the Implantation Sheath (IS) and the implant with fluoroscopy.
- Slowly advance the pusher (P) to unfold the distal disc in the aorta. Fig. 15.
- By pulling the Implantation Sheath (IS), gently pull the retention disc against the aortic ampulla. Fig. 16.


Fig. 16


Fig. 17

- Keep tension on the pusher cable while slowly withdrawing the Implantation Sheath (IS). At this point you can see the final configuration of the Nit-Occlud® PDA-R (D). Fig. 17.
- Verify the position of the device using a contrast injection through a pigtail catheter in the Aorta.
- Cut or remove the security seal (SS). Fig. 18.


Fig. 18


Fig. 19

- Hold firmly with one hand the “Y” connector and with the other pull the torquer at least 30 cm until the release of the device can be seen. Fig. 19 - 21.


Fig. 20


Fig. 21

- Perform a final angiographic evaluation 10 minutes after releasing the device.

Troubleshooting

- In case that the device cannot be released by pulling the torquer (T), it is suggested to elongate the pusher 2 cm near the “Y” connector. This will loosen the retaining wires and the device (D) will be released. Fig. 22 - 25.


Fig. 22: Position 1A


Fig. 24: Position 1B


Fig. 23: Position 2A


Fig. 25: Position 2B

Advice to Patients

- It is important that the patient reports possible allergies, especially to contrast mediums.
- It is recommended to give verbal and written information concerning the intervention and its possible complications.
- Bed rests 24 hours post implantation and relative rest during the next 90 days is recommended. Avoid competitive sports for two weeks.
- Avoid contacts with other individuals who could represent an infections risk.
- The clinical and echocardiographic follow-up of the patient is recommended at 24 hours, 30 days, 6 months and at 1 year.

Disposal after Use

Used medical products and accessories pose a potential biological hazard. For this reason, the products and their accessories should be handled and disposed of in accordance with recognised medical procedure and in compliance with the relevant legal regulations and local ordinances.

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